

1. Blanket Warmers

- Basis of Design
  - Steris (in-patient areas)
  - Blickman and Petigo (off-site outpatient clinics)
- Considerations
  - Full size units shall be on wheels to aid cleaning.
  - Half-size units must have a substantial base to sit upon
    - i. Unit needs to be on a stable surface.
    - ii. Unit is not to be larger than the surface it sets on.
    - iii. Need to be able to clean under the unit.
    - iv. Make sure the surface is strong to hold the weight of the warmer and blankets.
- Confirm sufficient electrical power is available. Minimize any other equipment connected to the circuit.

END OF SECTION